ULTRASTOP 15 ml Spray

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INSTRUCTIONS FOR USE

ULTRASTOP 15 ml Spray

ULTRASTOP 25 ml bottle with dropper pipette

ANTIFOG AGENT

INGREDIENTS / PROPERTIES

ULTRASTOP contains a special combination of surface active components in an alcohol/water matrix (ethanol, purified water, < 5 % surface-active substances).

ULTRASTOP is a high-quality anti-fogging agent for medical lenses. Examination mirrors and optical devices, as for example microscope oculars and lenses can be easily, simply and securely protected from fogging.

ULTRASTOP can also be used as anti-fogging agent in non-medical areas. All types of optics, goggles, glasses, lenses, mirrors etc. get quickly and reliably protected.

The special property of **ULTRASTOP** is that it reliably prevents fogging on the surfaces of medical lenses. Surface-active substances in a special combination form an absolutely clear and transparent film on the treated surfaces. This extremely thin film prevents the formation of droplets when water vapour becomes condensed; thus it prevents fogging of the surface.

AREAS OF USE

♦ For transient use (normally for less than 60 minutes) for all types of medical examination mirrors and lenses (e.g. larynx and epipharynx mirrors, colonoscopes, dental mirrors, etc.). Pre-warming the optics is no longer necessary as they can be pretreated with **ULTRASTOP** prior to consulting hours or home visits, which is time- and labour-saving for the doctor.

♦ For all types of glasses (e.g. operation glasses, work glasses, sun glasses, sports glasses, etc.).

♦ For all types of lenses and mirrors in microscopes or other scientific devices.

♦ For all other glass or mirror surfaces (diving masks, crash helmet visors, car glasses, etc.).

For use in strictly **aseptic endoscopy** the sterile form **ULTRASTOP "pro med."** (30 ml vial or 5 ml Monodose) has to be used.

Note: **ULTRASTOP** can be used on glass and synthetic materials. If other coatings have been used or if the coating is of doubtful resistance, please contact the manufacturer or the retailer of the object to be treated.

INSTRUCTIONS FOR USE

Apply **ULTRASTOP** directly or by means of an **ULTRASTOP**-moistened cloth to the dry and clean surface and spread evenly. Allow to take effect for a short time and after evaporation, if necessary, gently polish with a dry, clean and non-threading cloth.

ULTRASTOP if used properly, does not smear. However, if too much of it is applied, the surface may appear cloudy. In this case the surface should be cleaned with water and the above mentioned procedure should be repeated.

⚠ Highly flammable!

Keep out of the reach of children!

STORAGE CONDITIONS: Do not store above 25 °C. Protect from light.

OBTAINABLE AS

ULTRASTOP Spray, 15 ml

ULTRASTOP bottle with dropper pipette (bulb contains natural rubber), 25 ml

ULTRASTOP "pro med." Solution, 30 ml vial (**sterile**)

ULTRASTOP pro med. Solution Monodose (**sterile**):

1 application kit or 20 application kits consisting of **ULTRASTOP pro med.** Solution Monodose, 1 x 5 ml + 1 Sugi® Endopad, sterile endopad;

30 x 5 ml **ULTRASTOP pro med.** Solution Monodose (Clinic)

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Produttore: MoNo chem-pharm Produkte GmbH, Leystraße 129, A-1200 Vienna
Distribuzione: SIGMAPHARM Arzneimittel GmbH, A-1200 Vienna