BSS PLUS™

STERILE INTRAOCULAR IRRIGATING SOLUTION
(A Physiological Salt Solution Enriched with Bicarbonate,
Dextrose and Glutathione)

ENGLISH.....	2-4
DANSK	5-7
NEDERLANDS	8-10
EESTI	11-13
SUOMI.....	14-16
FRENCH.....	17-19
DEUTSCH	20-22
GREEK	23-26
MAGYAR	27-29
ITALIAN	30-33
LATVIEŠU	34-36
LIETUVIŲ	37-39
NORSK.....	40-42
POLSKI.....	43-46
ROMÂNĂ.....	47-49
PORTUGUÊS.....	50-52
ESPAÑOL	53-56
SVENSKA.....	57-59
УКРАЇНСЬКА.....	60-63
TURKISH.....	64-66

BSS PLUS™**STERILE INTRAOCULAR IRRIGATING SOLUTION**

(A Physiological Salt Solution Enriched with Bicarbonate, Dextrose and Glutathione)

DESCRIPTION

BSS PLUS™ Sterile Irrigating Solution is a weakly-buffered, sterile, physiologically balanced salt solution for use during all intraocular surgical procedures including those requiring a relatively long intraocular perfusion time. The electrolyte composition, pH, and osmolality approximate that of human aqueous humor. It does not contain a preservative and consists of two separate parts, which are aseptically mixed using the provided irrigating solutions transfer device just prior to use in surgery. The components of **BSS PLUS™** solution are biocompatible with the eye and have no pharmacological action.

COMPONENTS AND COMPOSITION PRIOR TO**RECONSTITUTION**

Part I: Part I is a sterile 480 mL solution in a 500 mL single-dose bottle to which the Part II concentrate is added. Each mL of Part I contains: Sodium Chloride 7.44 mg, Potassium Chloride 0.395 mg, Dibasic Sodium Phosphate 0.433 mg, Sodium Bicarbonate 2.19 mg, Hydrochloric Acid and/or Sodium Hydroxide (to adjust pH), in Water for Injection.

Part II: Part II is a sterile concentrate in a 20 mL single-dose vial for addition to Part I. Each mL of Part II contains: Calcium Chloride Dihydrate 3.85 mg, Magnesium Chloride Hexahydrate 5 mg, Dextrose 23 mg, Glutathione Disulfide (Oxidized Glutathione) 4.6 mg, in Water for Injection.

Mini Spike Transfer Device: The mini spike transfer device allows the transfer of **BSS PLUS™** Solution Part II into the glass bottle containing **BSS PLUS™** Solution Part I. The sterilized spike is made of polyethylene plastic and contains a 5 micron nylon filter inside the fluid path. The attached needle is 19 gauge x 1 inch B-Bevel. The transfer device is manufactured by B. Braun, and information for B. Braun is provided on the carton.

RECONSTITUTED PRODUCT COMPOSITION

After addition of **BSS PLUS™** Solution Part II to the Part I bottle, each mL of the reconstituted product contains: Sodium Chloride 7.14 mg, Potassium Chloride 0.38 mg, Calcium Chloride Dihydrate 0.154 mg, Magnesium Chloride Hexahydrate 0.2 mg, Dibasic Sodium Phosphate 0.42 mg, Sodium Bicarbonate 2.1 mg, Dextrose 0.92 mg, Glutathione Disulfide (Oxidized Glutathione) 0.184 mg, Hydrochloric Acid and/or Sodium Hydroxide (to adjust pH), in Water for Injection. The reconstituted product has a pH of approximately 7.4. Osmolality is approximately 305 mOsm. **BSS PLUS™** Solution is a moderately buffered solution containing a bicarbonate/phosphate buffering system. Acid Buffer Capacity = 0.015 M/L/pH; Base Buffer Capacity = 0.0024 M/L/pH.

INTENDED USE

BSS PLUS™ Sterile Irrigating Solution is intended for use by a trained ophthalmic surgeon in patients requiring intraocular surgical procedures, including those requiring a relatively long intraocular perfusion time (e.g. pars plana vitrectomy, phacoemulsification, extracapsular cataract extraction/lens aspiration, anterior segment reconstruction).

INDICATIONS

BSS PLUS™ Sterile Irrigating Solution is indicated for use as an intraocular irrigating solution during intraocular surgical procedures involving perfusion of the eye.

CONTRAINDICATIONS

There are no specific contraindications to the use of **BSS PLUS™** Solution when used as recommended. However, contraindications for the surgical procedure during which **BSS PLUS™** Solution is to be used should be strictly adhered to.

WARNINGS

1. This product is for IRRIGATION during ophthalmic surgery only. DO NOT use for injection or intravenous infusion.
2. The safety and effectiveness of using additives other than **BSS PLUS™** Part II have not been substantiated by Alcon.
3. DO NOT use unless product is clear, seal is intact, vacuum is present (Part I solution), and container is undamaged.
4. DO NOT use if sterile barrier is compromised or unintentionally opened before reconstitution.
5. DO NOT use if product is discolored or contains a precipitate.
6. DO NOT re-use this product. This product is for SINGLE patient use only and should not be used in more than one patient. Reuse of this single-use device may result in serious injury, such as but not limited to infection.
7. This solution contains no preservative; **DISCARD ANY UNUSED PORTION SIX HOURS AFTER PREPARATION.**

PRECAUTIONS

1. **DO NOT USE BSS PLUS™ SOLUTION UNTIL PART I IS FULLY RECONSTITUTED WITH PART II.**
2. Studies suggest that intraocular irrigating solutions which are iso-osmotic with normal aqueous fluids should be used with caution in diabetic patients undergoing vitrectomy since intraoperative lens changes have been observed.
3. There have been reports of corneal clouding or edema following ocular surgery in which **BSS PLUS™** was used as an irrigating solution. As in all surgical procedures appropriate measures should be taken to minimize trauma to the cornea and other ocular tissues.

PREPARATION

Reconstitute BSS PLUS™ Sterile Intraocular Irrigating Solution just prior to use in surgery by following the RECONSTITUTION DIRECTIONS provided. Strict aseptic procedures should be used during the reconstitution of BSS PLUS™ solution with the mini spike transfer device provided.

CLINICAL APPLICATIONS

None of the components of **BSS PLUS™** solution are foreign to the eye, and **BSS PLUS™** solution has no pharmacological action. Human perfused cornea studies have shown **BSS PLUS™** solution to be an effective irrigation solution for providing corneal detumescence and maintaining corneal endothelial integrity during intraocular perfusion. An *in vivo* study in rabbits has shown that **BSS PLUS™** solution is more suitable than normal saline or physiological salt solution for intravitreal irrigation because **BSS PLUS™** solution contains the appropriate bicarbonate, pH, and ionic composition necessary for the maintenance of normal retinal electrical activity. Human *in vivo* studies have demonstrated **BSS PLUS™** solution to be safe and effective when used during surgical procedures such as pars plana vitrectomy, phacoemulsification, cataract extraction/lens aspiration, and anterior segment reconstruction. No differences have been observed between adults and pediatric patients following use of this product.

GERIATRIC USE

No overall differences in safety or effectiveness have been observed between elderly and younger patients.

PERFORMANCE CHARACTERISTICS

A clinical evaluation of published market experience data, clinical investigations, and clinical literature data support that BSS PLUS™ Sterile Irrigating Solution has an acceptable safety and performance profile when acting as an intraocular irrigation solution during ophthalmic surgery. Adverse events with BSS PLUS™ Sterile Irrigating Solution occurred at a low rate, were mostly nonserious, and were expected when using an irrigating solution during intraocular surgery. The clinical benefits of using BSS PLUS™ Sterile Irrigating Solution for irrigation during ophthalmic surgery are greater than the risks associated with its use as BSS PLUS™ Solution is safe and performs as intended when used as directed.

CLINICAL BENEFITS

Based on a clinical evaluation of BSS PLUS™ Solution, the clinical benefits of BSS PLUS™ Sterile Intraocular Irrigating Solution in intraocular surgery include:

- Moistening of ocular tissues
- Protection of the corneal endothelium and other ocular tissues during intraocular surgeries
- Aid in maintaining intraocular spaces
- Rinsing away surgical debris
- Provides a suitable environment for the maintenance of normal structure and function of ocular tissues during prolonged ocular surgery
- Minimize potential complications during intraocular surgeries

ADVERSE REACTIONS

Postoperative inflammatory reactions, as well as incidents of corneal edema, increased intraocular pressure, and corneal decompensation, have been reported. Their relationship to the use of BSS PLUS™ Solution has not been established.

SERIOUS INCIDENT REPORTING

Any serious incidents that may reasonably be regarded as device related should be reported to Alcon Laboratories, Inc.:

By Phone:

Contact the local country office or your Alcon distributor.

By Email:

qa.complaints@alcon.com

Website:

<https://www.alcon.com/contact-us/>

Each device is identified by batch code which provides traceability, and this information should be given to Alcon.

NOTE: These serious incidents must also be reported to the competent authority for medical devices of the appropriate State.

OVERDOSAGE

The solution has no pharmacological action and thus no potential for overdosage. However, as with any intraocular surgical procedure, the duration of intraocular manipulation should be kept to a minimum.

DOSEAGE AND ADMINISTRATION

The solution should be used according to the standard technique employed by the operating surgeon. Use an administration set with an air-inlet in the plastic spike since the bottle does not contain a separate airway tube. Follow the directions for the particular administration set to be used. Insert the spike aseptically into the bottle through the target area of the rubber stopper. Allow the fluid to flow to remove air from the tubing before intraocular irrigation begins. If a second bottle is necessary to complete the surgical procedure, ensure that the vacuum is vented from the second bottle BEFORE attachment to the administration set. No special disposal of the BSS PLUS™ Sterile Irrigating Solution is required after use.

HOW SUPPLIED

BSS PLUS™ solution is supplied in two packages for reconstitution prior to use: a 500 mL bottle containing 480 mL (Part I) sterilized by steam and a 20 mL vial (Part II) sterilized by filtration. In addition, a mini spike transfer device sterilized by ethylene oxide is provided.

See the **PRECAUTIONS** section regarding reconstitution of the solution. Sterility is guaranteed until the expiration date unless the sterile barrier is damaged or opened.

STORAGE

Store Part I, Part II, and the mini spike transfer device between 2°C to 25°C. Discard prepared solution after six hours.

Alcon Laboratories, Inc.
6201 South Freeway
Fort Worth, Texas USA
76134-2099

EC REP

Alcon Laboratories Belgium
Lichterveld 3
2870 Puurs-Sint-Amands
Belgium

BSS PLUS™ STERILE INTRAOCULAR IRRIGATING SOLUTION

(A Physiological Salt Solution Enriched with Bicarbonate, Dextrose and Glutathione)

RECONSTITUTION DIRECTIONS: Use Strict Aseptic Technique

1. Remove the blue flip-off seal from the BSS PLUS™ Part I (480 mL) bottle. Remove the blue flip-off seal from the BSS PLUS™ Part II (20 mL) vial. Prepare the stoppers on both parts by using sterile alcohol wipes.
2. Peel open a BSS PLUS™ solution mini spike transfer device package (supplied), and remove the sterile transfer spike.
- NOTE:** This device is vented, permitting air to enter vial during solution transfer, thereby preventing the creation of a vacuum inside the vial. An air-inlet filter is provided to protect the system. Do not remove the air-inlet filter.
3. Remove protector from the white plastic piercing pin.
4. Firmly grasp device from behind the flange and insert the white plastic piercing pin into the upright rubber stopper of the BSS PLUS™ Part II (20 mL) vial.
5. Remove guard from filter needle. Firmly grasp vial in the palm of one hand, and with thumb and index finger, hold plastic flange against top of vial.
6. Invert vial and immediately insert filter needle into the rubber injection site of the BSS PLUS™ Part I (480 mL) bottle. (See illustration.)
7. Fluid will automatically transfer from the vial into the large vacuum bottle unless filter becomes occluded or loss of vacuum occurs.

NOTE: An excess amount of BSS PLUS™ Part II is provided in each vial. A non-transferred solution residual can be expected to remain in the vial.


8. Immediately remove needle from the BSS PLUS™ Part I container, and discard the needle after solution transfer has been completed (see Step 10).
 9. Mix the solution gently until uniform. Peel off the right-hand side of Part I bottle label (fully reconstituted BSS PLUS™ Solution). Record the patient's name **¶**? and the date **¶** and time **¶** of reconstitution. BSS PLUS™ Solution is now ready for use.
 10. Destroy needle from the mini spike transfer device after use. Use the guard to bend and break needle, and then, properly dispose of the used device in a sharps container.
- NOTE:** If the solution is not going to be used immediately, maintain aseptic conditions, and wipe stopper by using sterile alcohol wipe prior to use.

CAUTION: Reconstituted BSS PLUS™ Solution must be used within six hours of mixing. Discard any solution which has aged beyond that time. Never use the same bottle of BSS PLUS™ Solution on more than one patient.

Alternative Transfer Method

If preferred, the contents of the BSS PLUS™ Part II component may be aspirated aseptically with an 18-gauge cannula attached to a syringe and then transferred into the Part I bottle.

**HANGER DEVICE MUST BE
ASSEMBLED PROPERLY PRIOR
TO USE.**

SYMBOLS USED ON LABELING

SYMBOL	ENGLISH
	Medical device
	Caution
	Consult instructions for use or consult electronic instructions for use
	Do not re-use
	Do not re-sterilize
	Temperature Limit
	Non-pyrogenic
	Sterilized using aseptic processing techniques
	Sterilized using steam

	Sterilized using ethylene oxide
	Single sterile barrier system
	Use-by date
	Manufacturer
	Date of manufacture
	Batch code
	Fragile, handle with care
	Authorised Representative in the European Community
	Distributor
	Not for intravenous use
	Do not use if package is damaged and consult instructions for use
	Patient identification
	Date
	Time
	Open here
	Country of Manufacture

Alcon Laboratories, Inc.
6201 South Freeway
Fort Worth, Texas
76134-2099 USA

Distributed in Australia by:
Alcon Laboratories
(Australia) Pty Ltd
15 Talavera Road
Macquarie Park NSW 2113
Australia
Phone Toll Free: 1800 224 153

EC REP

Alcon Laboratories Belgium
Lichterveld 3
2870 Puurs-Sint-Amands
Belgium

Distributed in New Zealand by:
Alcon Laboratories
(New Zealand) Limited
12 St Marks Road
Remuera, Auckland 1050
New Zealand
Free Phone: 0800 101 106

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
HANGER DEVICE MUST BE
ASSEMBLED PROPERLY PRIOR
TO USE.

SYMBOLS USED ON LABELING

SYMBOL	ENGLISH
	Medical device
	Caution
	Consult instructions for use or consult electronic instructions for use
	Do not re-use
	Do not resterilize

	Temperature Limit
	Non-pyrogenic
	Sterilized using aseptic processing techniques
	Sterilized using steam
	Sterilized using ethylene oxide
	Single sterile barrier system
	Use-by date
	Manufacturer
	Date of manufacture
	Batch code
	Fragile, handle with care
	Authorised Representative in the European Community
	Distributor
	Not for intravenous use
	Do not use if package is damaged and consult instructions for use
	Patient identification
	Date
	Time
	Open here
	Country of Manufacture

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