Codman

**CODMAN® Refill Kit (AS 10)
for ARCHIMEDES® Pumps**

EC REP

Codman

A division of Johnson & Johnson Medical Ltd.
Pinewood Campus, Nine Mile Ride
Wokingham RG40 3EW
United Kingdom

CE
0086

Authorized 2003

LCN 200370-001/K

© 2003-2011 Codman & Shurtleff, Inc.

 Revised 10/11


Table of Contents

Indications	4
Contraindications	5
Warnings	6
Precautions	7
Adverse Events	8
Drug Stability Information	9
Product Description	10
Sterility	11
Instructions for Use	12
Refilling the Implanted Pump	12
A. Setting Up and Priming the Refill Assembly	13
B. Emptying the Drug Reservoir of the Remaining Drug	14
C. Rinsing the Drug Reservoir	15
D. Refilling the Drug Reservoir	16
E. Completing the Documentation	18
Scheduling a Pump Refill	18
Discontinuing Drug Therapy	19
Warranty	20
Appendix A: Emergency Procedures	21
Drug Overdose	21
Appendix B: Calculating the Drug Infusion Solution	22
Appendix C: Troubleshooting	24
Troubleshooting Table	24
Troubleshooting Flowchart	27

Indications

ENGLISH

IMPORTANT INFORMATION

Please Read Before Use

CODMAN® Refill Kit (AS 10) for ARCHIMEDES® Pumps

STERILE EO

Become thoroughly familiar with the information contained in this user manual prior to filling or refilling an ARCHIMEDES Implantable Infusion Pump. Failure to follow these instructions can result in patient complications ranging from a failure of the intended therapy to a drug under- or overdose.

For additional information, refer to the user manual packaged with the infusion pump.

Indications

The CODMAN Refill Kit (AS 10) is used only to fill, refill, or empty the drug reservoir of an ARCHIMEDES Implantable Infusion Pump. The pumps are indicated for the intrathecal infusion of the drugs listed below.

- For use with preservative-free morphine sulfate sterile solution (4.0–6.5 pH) in the treatment of chronic pain
- For use with LIORESAL® Intrathecal (baclofen) (5.0–7.0 pH) for relief of spasticity

Use a 0.9% solution of preservative-free sodium chloride injection (4.5–7.0 pH) to achieve the physician-prescribed drug concentration.

WARNING: Use only the drugs and diluent listed above with the ARCHIMEDES pump.

Using incompatible or contaminated (out of specification) drug solutions can damage the pump flow restrictor chip and cause an increase in flow rate. This can lead to serious complications, including overdose and death.

WARNING: Refer to the appropriate drug labeling for a complete list of drug indications, contraindications, warnings, precautions, adverse reactions, dosage and administration information, screening procedures, and overdose procedures. Physicians prescribing an ARCHIMEDES Implantable Infusion Pump for use with these drugs must be familiar with the drug stability information listed in *Drug Stability Information*.

WARNING: Based on studies carried out by Codman that evaluated precipitation of Morphine Sulfate for Intrathecal Injection at varying pH levels, the use of preservative-free morphine sulfate solutions with pH levels above 6.5 is not recommended for these pumps. In particular, we have found that higher concentrations of morphine sulfate solutions (>10 mg/mL) might be more susceptible to precipitation when the pH levels of the drug are at or near the upper pH limits specified in USP requirements for this drug (i.e., at or near pH 6.5).

Precipitation of morphine could result in an accumulation of the precipitates in the flow path causing possible slowing or stoppage of flow through the flow-restrictor chip. This can lead to inadequate pain relief and/or the need to remove and replace the pump.

Contraindications

There are no known contraindications for the use of the CODMAN Refill Kit (AS 10).

Observe all contraindications relating to the use of the prescribed drug.

Warnings

WARNINGS

Improper use of implanted infusion pumps can result in drug under- or overdose. Users must comply with product instructions for initial filling and refilling of the pump. **Only qualified medical personnel must perform these procedures.**

Do not mix drugs. The effects of mixing drugs in the pump are not known.

Use only the drugs and diluent listed in *Indications* with the ARCHIMEDES pump. Using incompatible or contaminated (out of specification) drug solutions can damage the pump flow restrictor chip and cause an increase in flow rate. This can lead to serious complications, including overdose and death.

The physician must consider the drug concentration, dose, and flow rate relationships of the pump before selecting the pump volume and flow rate. Failure to consider these factors can lead to a drug under- or overdose.

Do not overfill the drug reservoir. Always identify the drug reservoir capacity (volume) before performing a filling procedure. Always follow the instructions for removing the remaining drug solution before refilling the drug reservoir. Failure to remove the remaining drug solution can result in overfilling the drug reservoir. Overfilling can result in overinfusion, which can lead to a drug overdose or cause damage to the pump.

It is very important to correctly calculate the infusion solution. Calculation errors can result in overinfusion or underinfusion. Refer to *Appendix B: Calculating the Drug Infusion Solution*.

It is very important to access the central port during refill procedures. An error in accessing the central port can result in tissue damage or a drug under- or overdose, or drug withdrawal symptoms.

Never inject into the pump catheter or the intraspinal catheter; a drug under- or overdose or damage to the catheter can result.

Do not use the needle provided with the refill kit for accessing the bolus port.

Injection in the bolus port or into the pump pocket can result in drug overdose. To prevent injection errors:

- identify the location of the central port and the bolus port;
- use the SURECAN® Filling Needle provided in the CODMAN Refill Kit (AS 10) for accessing the central port;
- confirm drug reflux into the syringe barrel during the filling procedure.

After rinsing the drug reservoir, intrathecal delivery of the new drug solution from the drug reservoir will be delayed until the solution traverses the total catheter length. Carefully calculate the amount of the delay, based on the total catheter length and the flow rate, before administering additional medication by another delivery method.

Warnings – 6

Precautions

Use sterile technique in all phases of handling this product.

Inspect the sterile package carefully. Do not use if:

- the package or seal appears damaged;
- the contents appear damaged; or
- the expiry date has passed.

The refill kit is for **single use only. Do not reuse.**

If local or systemic infection is suspected, use extreme caution when emptying and/or refilling the drug reservoir. If infection is suspected perform appropriate diagnostic procedures and intervention.

Use only the SURECAN Filling Needle provided with the CODMAN Refill Kit (AS 10) for filling and refilling ARCHIMEDES Infusion Pumps. This needle is designed specifically for use with these pumps. Use of other needles can damage the central port and/or result in a failure to administer solution appropriately.

Use only preservative-free solutions for intrathecal applications.

Adverse Events

Adverse Events

Adverse events related to emptying, filling, or refilling the drug reservoir include:

- Contamination of the drug reservoir resulting in infection, such as meningitis;
- Injection errors, which can lead to tissue damage or a drug under- or overdose; or
- Overfilling of the drug reservoir, which can lead to drug overdose or damage to the pump.

Complications relating to the drug therapy, such as:

- drug toxicity (local and systemic) and its related side effects;
- complications due to the use of a drug that is not approved for intrathecal administration;
- complications due to the use of drugs that are not approved for use with the system;
- complications due to use of an approved drug or diluent that is out of specification;
- extravasation.

Drug Stability Information

Use only the drugs listed below with ARCHIMEDES Infusion Pumps for the indications described in *Indications*. Mix these drugs with only 0.9% solution of preservative-free sodium chloride injection (sterile 0.9% saline solution) to achieve the appropriate concentration.

Drug	Stability
Preservative-free morphine sulfate	120 days
LIORESAL Intrathecal (baclofen)	120 days

Refer to the drug labeling for complete prescribing information, indications, contraindications, warnings, precautions, and adverse events.


When filling the drug reservoir, verify that the “Use By” date of the drug will not occur before the patient’s next refill.

Product Description

Product Description

The CODMAN Refill Kit (AS 10) (Figure 1) contains the following components, provided sterile:

- 2 SURECAN Filling Needles with grey luer-lock connector
- 10-mL collection syringe without plunger, silicone stopper included
- 2 connecting tubes
- 3-way stopcock
- 60-mL syringe
- 10-mL syringe
- Filter, 0.2 micron
- Fenestrated drape (not shown)

Product Description – 10

Sterility

②

This product is for SINGLE USE ONLY; DO NOT RESTERILIZE.

Codman Single Use devices have not been designed to undergo or withstand any form of alteration, such as disassembly, cleaning or resterilization, after a single patient use. Reuse can potentially compromise device performance and any usage beyond the design intent of this single-use device may result in unpredictable loss of functionality.

Codman & Shurtleff, Inc. will not be responsible for any product that is resterilized, nor accept for credit or exchange any product that has been opened but not used.

As long as the inner unit is not opened or damaged, the product is sterile.

The following components have been tested and were determined to be nonpyrogenic:

SURECAN Filling Needles
Luer-lock connector
Connecting tubes

The following components have been tested and were determined to have a nonpyrogenic fluid pathway:

10-mL collection syringe
60-mL syringe
10-mL syringe
Filter
3-way stopcock

Instructions for Use – Refilling the Implanted Pump

INSTRUCTIONS FOR USE

These instructions do not address the procedure for initial filling of the drug reservoir in preparation for implantation. Follow the instructions in *Preparing the Pump* in the user manual that is packaged with the ARCHIMEDES Implantable Infusion Pump.

Refilling the Implanted Pump

PRECAUTIONS

Ensure that all luer-lock connections are secure.

Never overfill the drug reservoir. Always be sure to identify the reservoir capacity before performing a refilling procedure.

Use only the SURECAN Filling Needle provided with the CODMAN Refill Kit (AS 10), to access the central port. Use of other needles can damage the central port and/or result in a failure to administer solution appropriately.

Always follow the instructions for removing the remaining drug solution before refilling the drug reservoir. Failure to remove the remaining drug solution can result in overfilling the reservoir. Overfilling can result in overinfusion, which can lead to a drug overdose or cause damage to the pump.

Do not aspirate during emptying of the drug reservoir. Allow pump pressure to evacuate the drug reservoir during the emptying procedure. Aspiration of the reservoir can damage the pump's drive system.

The drug reservoir contents are under significant pressure. To prevent the reservoir contents from being ejected, do not use an open syringe when emptying the drug reservoir.

Always ensure proper needle placement (needle held perpendicular to the pump and inserted completely to the needle stop) before refilling the reservoir.

Use of excessive force when inserting the filling needle into the central port can damage the needle tip. Never use a filling needle if the tip is bent; use the second filling needle provided. Using a bent filling needle will cause damage to the central port, and can result in a failure to administer solution appropriately.

The procedure for refilling an implanted pump consists of five steps:

- A. Setting up and priming the refill assembly;
- B. Emptying the drug reservoir of the remaining drug;
- C. Rinsing the drug reservoir (required to change the drug or drug concentration);
- D. Refilling the drug reservoir; and
- E. Completing the documentation.

Required materials:

- Aspiration needle
- CODMAN Refill Kit (AS 10)
- Appropriate amount of drug solution (volume to be determined by the capacity of the drug reservoir)
- 3 mL of drug solution to prime the refill assembly
- Skin disinfectant
- Self-adhesive sterile bandage (optional)

Additional materials if rinsing the drug reservoir:


- Aspiration needle
- (2) 10-mL syringes with luer-lock connectors
- 10 mL of sterile 0.9% saline solution

A. Setting Up and Priming the Refill Assembly (see Figures 2a, 2b and 2c)


1. Check the patient's records for the drug reservoir capacity and verify that the amount of drug solution does not exceed the drug reservoir capacity plus the approximately 3 mL lost in priming.
2. Use the aspiration needle to fill the 60-mL syringe with the drug solution. Expel the air from the syringe. **Remove the aspiration needle.**
3. Assemble the parts (see Figure 2a).
 - a. Connect the filter to the 60-mL syringe.
 - b. Connect one of the connecting tubes to the filter.
 - c. Connect the free end of this connecting tube to the **perpendicular fitting of the stopcock.**
 - d. Connect the SURECAN Filling Needle to the second connecting tube.
 - e. Connect the second connecting tube to the stopcock male fitting.
 - f. Turn the stopcock handle so it points to the filling needle.
 - g. Slowly depress the plunger of the 60-mL syringe until drug solution appears at the open female fitting of the stopcock.
 - h. Ensure that the plunger of the 10-mL syringe is fully depressed. Connect the 10-mL syringe to the stopcock.
 - i. Fill the 10-mL syringe with 5 mL of solution by slowly depressing the plunger of the 60-mL syringe (see Figure 2b).
 - j. Turn the stopcock handle so it points to the perpendicular fitting. If air is present in the 10-mL syringe, remove the syringe and expel the air; reconnect the syringe to the stopcock.
 - k. Remove the protective cover from the filling needle.
 - l. **Slowly** depress the plunger of the 10-mL syringe until drug solution exits from the filling needle (Figure 2c). The assembly is now primed. Replace the protective cover on the filling needle.
4. Turn the stopcock handle so it points to the filling needle. Remove the 10-mL syringe. Place the 10-mL syringe and the remainder of the assembly on a sterile field.

Note: Once the system is primed, monitor the dosage amount through the use of the 10-mL syringe. Deliver only the amount that corresponds to the drug reservoir capacity. If less than 3 mL was used to prime the system, discard the excess in the appropriate manner.


2a

2b

2c

Instructions for Use – Refilling the Implanted Pump

B. Emptying the Drug Reservoir of the Remaining Drug (see Figure 3)

CAUTION: Always ensure that the stopcock handle is positioned correctly, as outlined in the procedure that follows. Positioning the stopcock handle incorrectly can result in an unintended effect, such as drug solution refluxing from the pump into the 60-mL syringe.

Refluxed fluid from the pump is of particular concern in the following instances:

- When changing the drug – because the new drug could become altered with refluxed drug
- When changing the concentration of the drug – because the new drug concentration could become altered when refluxed drug solution is mixed with the new

If either of these instances occurs, discard the refill kit and drug solution in the appropriate manner. Use a new refill kit and new quantity of drug solution.

Perform the following calculation to determine the volume of drug solution removed from the pump, if the pump contents reflux into the 60-mL syringe:

$$\begin{aligned} & \text{Volume (mL) of drug solution in the 60-mL syringe} \\ & + \text{Volume (mL) of drug solution in the 10-mL syringe} \\ & + 3 \text{ mL (amount of drug solution used to prime the system)} \\ & - \text{Original amount (mL) of new drug solution} \\ & = \text{Amount of drug solution refluxed from the pump (mL)}^* \end{aligned}$$

*Use this amount in Step 2e and Step 3a, which follow.

1. Set up
 - a. Palpate the pump. Locate both the central port and the bolus port.
 - b. Disinfect a wide area of skin over and around the **central port**; cover the prepared area with the sterile fenestrated drape.
2. Empty the drug reservoir
 - a. Connect the **collection syringe** onto the stopcock. Firmly place the silicone stopper onto the end of the collection syringe barrel.
 - b. Ensure that the stopcock handle points to the filling needle. Remove the protective cover from the filling needle. Holding the filling needle perpendicular to the pump, puncture the central port. Insert the filling needle until it contacts the needle stop. **CAUTION: If the needle does not contact the needle stop, it will not be possible to remove drug from the pump.**

- c. Loosen the silicone stopper and hold the collection syringe upright. Turn the stopcock handle so it points to the perpendicular fitting. The contents of the reservoir will reflux into the collection syringe. If reflux does not appear, check the position of the filling needle. If needed, twist the needle to ensure that it penetrates the central port. **If, after checking the filling position, no reflux is present, consult Appendix C: Troubleshooting.**
 - d. When the reflux of drug has stopped, wait an additional 30 seconds to ensure that the pump is empty. **Turn the stopcock handle so it points to the filling needle.**
 - e. Record the volume removed from the drug reservoir.
 - f. Secure the silicone stopper firmly to prevent spillage. Turn the collection syringe upside down before disconnecting it from the connecting tube to prevent small amounts of spillage.
 - g. Remove the collection syringe from the stopcock. The filling needle remains in place in the central port; the tubing and remainder of the refill assembly are attached to the needle.
3. Calculate the flow rate
 - a. Calculate the flow rate according to the following equation:


$$\frac{\text{Refilling volume from previous visit (mL)} - \text{Volume removed from current visit (mL)}}{\text{Running time since last refill (days) [24 h]}} = \text{Calculated flow rate (mL/24 h)}$$

If the difference between the calculated flow rate and the set flow rate is unacceptable, refer to Appendix C: Troubleshooting.

- b. Dispose of the syringe and refluxed solution appropriately.
- c. If there will be a change in the drug solution (either the drug itself or the concentration), proceed to *Rinsing the Drug Reservoir*. If there will be no change to either the drug or the concentration, proceed to *Refilling the Drug Reservoir*.

Instructions for Use – 14

3

C. Rinsing the Drug Reservoir

WARNING: Use caution to prevent drug overdose when changing concentrations in the drug reservoir. Before changing the concentration of the drug solution, first empty the drug reservoir, then rinse the reservoir twice with sterile 0.9% saline solution as instructed.

This procedure is performed with the filling needle still in place in the central port.

1. Using an aspiration needle, draw 5 mL of saline into each of two 10-mL syringes.
2. Connect one saline-filled syringe to the straight female fitting of the stopcock. **Turn the stopcock handle so it points to the perpendicular fitting.**
3. Slowly inject the saline into the reservoir. Remove pressure from the syringe plunger and allow the saline to reflux into the syringe. Remove and discard the syringe and the refluxed solution appropriately.
4. Connect the second saline-filled syringe to the stopcock.
5. Slowly inject the saline into the reservoir. Remove pressure from the syringe plunger and allow the saline to reflux into the syringe. **Turn the stopcock handle so it points to the filling needle.**
6. Remove and discard the syringe and the refluxed solution appropriately.
7. Leave the filling needle in place in the central port.
8. Proceed to *Refilling the Drug Reservoir*.

Instructions for Use – Refilling the Implanted Pump

D. Refilling the Drug Reservoir (refer to Figures 4a and 4b)

CAUTION: Before refilling the reservoir, ensure that the remaining drug solution has been removed (see Emptying the Drug Reservoir).

CAUTION: Before changing the concentration of the drug solution, first empty the reservoir, then rinse the reservoir twice with sterile 0.9% saline solution as instructed in Rinsing the Drug Reservoir.

CAUTION: The contents of the pump are under pressure. In performing Steps 3b and 3c, do not completely remove your thumb from the plunger of the 10-mL syringe. Use your thumb to keep the plunger under control. Failure to control the plunger may cause drug solution to spill.


Note: Once the system is primed, monitor the dosage amount through the use of the 10-mL syringe. Deliver only the amount that corresponds to the drug reservoir capacity. If less than 3 mL was used to prime the system, discard the excess in the appropriate manner.

1. Ensure that the stopcock handle is turned so it points to the filling needle. Reattach the 10-mL syringe containing 5 mL of drug solution to the straight female fitting of the stopcock.
2. Check the drive function of the pump.
 - a. Depress the plunger of the 60-mL syringe until the 10-mL syringe fills with drug solution.
 - b. **While keeping your thumb on the plunger of the 10-mL syringe ready to counteract the reflux pressure of the pump**, turn the stopcock handle so it points to the perpendicular fitting. Inject the contents of the 10-mL syringe into the reservoir.
 - c. Check the function of the pump's drive system by releasing the pressure from the plunger of the 10-mL syringe for a few seconds. Reflux **must** appear in the syringe.
 - d. Inject the remaining contents of the 10-mL syringe into the reservoir.
 - e. **While continuing to exert pressure on the 10-mL syringe plunger, turn the stopcock handle so it points to the filling needle.**
3. Inject the remaining drug solution.
 - a. Depress the plunger of the 60-mL syringe until the 10-mL syringe fills with drug solution.
 - b. **While keeping your thumb on the plunger of the 10-mL syringe ready to counteract the reflux pressure of the pump**, turn the stopcock handle so it points to the perpendicular fitting. Inject the contents of the 10-mL syringe into the reservoir.
 - c. **While continuing to exert pressure on the 10-mL syringe plunger, turn the stopcock handle so it points to the filling needle.**
 - d. Repeat Steps 3a through 3c until the correct volume of the drug solution has been injected.
4. Remove the filling needle from the central port and, if necessary, cover the injection site with a self-adhesive sterile bandage.


Instructions for Use – 16

Instructions for Use – Refilling the Implanted Pump

4a

4b

Instructions for Use – Scheduling a Pump Refill

E. Completing the Documentation

Complete the patient passport.

Scheduling a Pump Refill

CAUTION: At refill, the pump must contain at least 3 mL of fluid. When scheduling the first refill, take into account the reservoir priming volume of approximately 1 mL. The flow rate of the pump decreases rapidly as the volume in the reservoir decreases below 3 mL. This can result in the potential loss of therapeutic effect or drug withdrawal symptoms.

Note: The flow rate is affected by changes in altitude and/or temperature. Refer to the diagrams in the user manual packaged with the ARCHIMEDES Pumps to determine the significance of the change.

1. Calculate the Refill Interval

In the equation below, the 3 mL subtracted from the reservoir capacity reflects the amount of drug solution that must remain in the drug reservoir to ensure that the flow rate remains constant.

$$\frac{\text{Reservoir capacity} - 3 \text{ mL}}{\text{Flow Rate}} = \text{Refill Interval}$$

Example:

Reservoir capacity: 20 mL

Flow rate: 0.5 mL/24 h

$$\frac{20 \text{ mL} - 3 \text{ mL}}{0.5 \text{ mL}/24 \text{ h}} = 34 \text{ days [24 h]}$$

In this example, the patient must be scheduled to return within 34 days.

2. Schedule the refill appointment with your patient.

Discontinuing Drug Therapy

CAUTION: If the drug therapy is to be discontinued for a prolonged period of time, empty the pump of the drug solution and refill it with sterile 0.9% saline solution, following the prescribed refill schedule to maintain a patent fluid pathway.

To discontinue drug therapy, follow the instructions for *Refilling the Implanted Pump* using saline in place of drug solution.

Warranty

WARRANTY

Codman & Shurtleff, Inc. warrants that this medical device is free from defects in both materials and workmanship. **Any other express or implied warranties, including warranties of merchantability or fitness, are hereby disclaimed. Suitability for use of this medical device for any particular surgical procedure should be determined by the user in conformance with the manufacturer's instructions for use. There are no warranties that extend beyond the description on the face hereof.**

® CODMAN and ARCHIMEDES are registered trademarks of Codman & Shurtleff, Inc.

® Lioresal is a registered trademark of Novartis Pharma AG

® SURECAN is a registered trademark of B. Braun Medical

APPENDIX A: EMERGENCY PROCEDURES

Drug Overdose

WARNING: REFER TO THE APPROPRIATE DRUG LABELING FOR A COMPLETE LIST OF DRUG INDICATIONS, CONTRAINDICATIONS, WARNINGS, PRECAUTIONS, ADVERSE REACTIONS, DOSAGE AND ADMINISTRATION INFORMATION, SCREENING PROCEDURES, AND OVERDOSE PROCEDURES CONTAINED IN THE PRESCRIBED DRUG LABELING.

Consult the patient's medical record or with the patient's physician to confirm the drug or drug concentration within the drug reservoir.

In the event of a suspected overdose:

1. Aspirate the catheter, if appropriate.
2. Empty the reservoir of the remaining drug. Refer to *Section B, Emptying the Drug Reservoir of the Remaining Drug*, to empty the pump reservoir.
3. After stabilizing the patient, rinse the reservoir. Refer to *Section C, Rinsing the Drug Reservoir*.
4. Once the drug reservoir is rinsed, fill the reservoir with **0.9%** saline solution to keep the pathway patent, or fill with the correct prescription.

Appendix B: Calculating the Drug Infusion Solution

APPENDIX B: CALCULATING THE DRUG INFUSION SOLUTION

The drugs indicated for use with the ARCHIMEDES Implantable Infusion Pump must be diluted only with sterile **0.9%** saline solution to make the infusion solution.

CAUTION: Preservative-free morphine sulfate is supplied in various concentrations. Be aware of the concentration before calculating the amount of drug to be infused.

In order to calculate the drug infusion solution, one must perform 3 steps:

- calculate the running time of the pump;
- calculate the amount of the drug required in [mg] or [μ g]; then convert into [mL]; and
- calculate the amount of saline required.

CAUTION: Prepare an additional 3 mL of drug solution, in excess of the capacity of the reservoir, to prime the refill assembly.

Calculate the running time of the pump:

$$\frac{\text{Reservoir capacity}}{\text{Flow Rate/24 hours}} = \text{Running time (days) [24 h]}$$

Calculate the amount of the drug:

$$\begin{array}{|c|} \hline \text{Required daily dosage} \\ \text{of the drug} \\ \hline \text{[mg/24 h] or } [\mu\text{g/24 h}] \\ \hline \end{array} \times \begin{array}{|c|} \hline \text{Running time of the} \\ \text{infusion pump in} \\ \hline \text{[24 h] days} \\ \hline \end{array} = \begin{array}{|c|} \hline \text{Total amount of} \\ \text{drug required} \\ \hline \text{[mg] or } [\mu\text{g}] \\ \hline \end{array}$$

Convert the total amount of drug required from mg or μ g to mL. In the case of drugs with lower dosages (such as baclofen), please observe the conversion factor for converting mg to μ g (1 mg = 1000 μ g).

$$\begin{array}{|c|} \hline \text{Total amount of drug} \\ \hline \text{[mg] or } [\mu\text{g}] \\ \hline \end{array} / \begin{array}{|c|} \hline \text{Concentration of drug} \\ \hline \text{[mg/mL] or } [\mu\text{g/mL}] \\ \hline \end{array} = \begin{array}{|c|} \hline \text{Total amount of} \\ \text{drug required} \\ \hline \text{[mL]} \\ \hline \end{array}$$

Calculate the amount of sterile 0.9% saline solution:

$$\begin{array}{|c|} \hline \text{Reservoir capacity} \\ \hline \text{[mL]} \\ \hline \end{array} - \begin{array}{|c|} \hline \text{Total amount of} \\ \text{drug required} \\ \hline \text{[mL]} \\ \hline \end{array} = \begin{array}{|c|} \hline \text{Total amount of} \\ \text{saline [mL]} \\ \hline \end{array}$$

Result: Total amount of infusion solution

$$\begin{array}{|c|} \hline \text{Total amount of} \\ \text{drug required} \\ \hline \text{[mL]} \\ \hline \end{array} + \begin{array}{|c|} \hline \text{Total amount of} \\ \text{saline [mL]} \\ \hline \end{array} = \begin{array}{|c|} \hline \text{Total amount of} \\ \text{infusion solution} \\ \hline \text{[mL]} \\ \hline \end{array}$$

Example of a Drug Calculation

Note: The preservative-free morphine sulfate concentration in this example is 10 mg/1 mL.

In this example:

- the daily dosage is 4 mg morphine
- the reservoir capacity is 20 mL
- the pump flow rate is 0.8 mL/24 hours

Step 1: Calculate the running time of the pump

$$\frac{20 \text{ mL reservoir capacity}}{0.8 \text{ mL flow rate/24 h}} = 25 \text{ days running time}$$

Step 2: Calculate the total amount of drug required

$$4\text{-mg daily dosage} \times 25 \text{ days running time} = 100 \text{ mg morphine}$$

Step 2a: Calculate the number of morphine ampules (assuming 1 mL per ampule and 10 mg morphine/mL)

$$\frac{100 \text{ mg morphine}}{10 \text{ mg/mL}} = 10 \text{ ampules @ 1 mL/ampule} = 10 \text{ mL morphine}$$

Step 3: Calculate the amount of saline

$$20 \text{ mL reservoir} - 10 \text{ mL morphine} = 10 \text{ mL saline}$$

Result

The infusion pump is filled with a mixture of:

10 mL	morphine mix (at 10 mg morphine/mL = 100 mg total morphine)
+ 10 mL	saline
= 20 mL	infusion solution

Appendix C: Troubleshooting

APPENDIX C: TROUBLESHOOTING

Troubleshooting Table

The following table lists symptoms that may occur, along with their possible causes, methods of verification, and suggested remedial actions. If, after following the suggested solution(s) you cannot restore pump function to normal, please contact your Codman representative.

Symptoms	Possible Cause	Method of Verification	Remedial Action
Patient symptoms reappear or intensify	• Pump leakage • Pump failure (no output) • Catheter occlusion • Catheter tear • Catheter flowpath interruption • Drug concentration incorrect	Refer to Troubleshooting Flowchart	Refer to Troubleshooting Flowchart
At refill, the amount of fluid in the reservoir exceeds the tolerances expected (more or less fluid than expected)	Pump failure	Check to see if patient symptoms have reappeared or intensified	Refer to Troubleshooting Flowchart
	Damage to the flow restrictor chip	Withdraw and replace the drug solution and evaluate the flow rate	Contact your Codman representative and consider elective replacement of the pump
	Extenuating circumstances	Examine the patient history to determine if the patient has experienced high or low altitudes, fevers, etc.	Empty the pump, refill, and monitor. Explain to patient the effects of high and low altitude, fever, etc. on pump performance
	Incorrect refill calculations	Recalculate the refill period	Refill pump as required, and give patient a new refill schedule
Unable to palpate pump ports	Pump turned upside down	Use radiography to verify pump position	Perform surgical revision to correct pump positioning
	Pump has rotated	Use radiography to verify pump position	Refill pump as required, and give patient a new refill schedule
	Swelling at pump site	Aspirate pump pocket or site swelling and verify contents	Treat according to aspiration content analysis
	Patient weight gain	Review patient's history	Inform patient of effects of weight gain on treatment

Appendix C: Troubleshooting

Swelling or tenderness at the pump site	Pump or catheter system leakage	Aspirate fluid at site and verify contents	If aspiration contents show pump leakage or catheter system leakage, perform the following: 1. Perform radiography to determine the source of the leak 2. Perform surgical revision on pump or catheter system as required
	Infection at site	Aspirate fluid at site and verify contents	Treat according to aspiration contents
	Foreign body rejection phenomena	Aspirate fluid at site and verify contents	Treat according to aspiration contents
	Twiddler's Syndrome		
Treatment performance is erratic	Disease related	Physician assessment	Treat accordingly
	Intermittent kinking of catheters	Use radiography to verify catheter integrity	Treat accordingly
	Extenuating circumstances	Examine the patient history to determine if the patient has experienced high or low altitudes, fevers, etc.	Explain to patient the effects of high and low altitude, fever, etc. on pump performance
Unable to inject bolus	Incorrect needle used	Check the needle to make sure the bolus needle provided is being used	Use only the bolus needle provided and retry
	The bolus needle is positioned incorrectly	Check the needle position	Hold needle perpendicular to the pump and retry
	The bolus needle is damaged or occluded	Inspect the bolus needle	If bolus needle is suspected to be damaged or occluded, discard and use a new bolus needle
	Pump is positioned incorrectly	Use radiography to check pump position and catheter system	Treat accordingly
	The catheter is bent or occluded	Use radiography to check pump position and catheter system	Treat accordingly

Appendix C: Troubleshooting

Unable to inject into central port	Incorrect needle used	Make sure the SURECAN Filling Needle is being used	Use only the SURECAN Filling Needle provided and retry
	The SURECAN Filling Needle is improperly positioned	Check the needle position	Ensure the needle is being held perpendicular to the central port and retry, twisting needle to ensure penetration
	The pump is positioned incorrectly	Use radiography to check pump position	Treat accordingly
	The SURECAN Filling Needle is damaged or occluded	Inspect the needle	If the needle appears damaged or occluded, discard and use a new needle

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