BRAUN

ExactFit™ 1

BUA 5000

GB English.....	2
DE Deutsch.....	14
ES Español.....	26
FR Français.....	38
HR Hrvatski.....	50
HU Magyar.....	62
PL Polski.....	74
RO Română.....	86
RU Русский.....	98
SI Slovenščina.....	110
TR Türkçe.....	122


www.hot-europe.com

This product is manufactured by Kaz Europe Sàrl under a license to the 'Braun' trademark. 'Braun' is a registered trademark of Braun GmbH, Kronberg, Germany.

 Kaz Europe Sàrl
Place Chauderon 18
CH-1003 Lausanne
Switzerland

© 2015, Kaz Europe Sàrl

Made and printed in China


BUA5000EU

Intended use of Braun ExactFit™ 1

Braun ExactFit™ 1 is a blood pressure monitor that has been developed for easy and accurate blood pressure measurements from the upper arm. The measuring accuracy of Braun ExactFit™ 1 was tested at the time of manufacture and was proven by clinical research in accordance with ESH.[†]

What you should know about blood pressure

Blood pressure constantly changes throughout the day. It rises sharply in the early morning and declines during the late morning. Blood pressure rises again in the afternoon and finally drops to a low level at night. Also, it may vary in a short period of time. Therefore, readings from successive measurements can fluctuate.

Note: Blood pressure readings from a healthy 31-year-old male, measured at 5-minute intervals

Blood pressure measured in a doctor's office only provides a momentary value. Repeated measurements at home better reflect one's actual blood pressure values under everyday conditions.

Moreover, many people have a different blood pressure when they measure at home, because they tend to be more relaxed than when in the doctor's office.

Regular blood pressure measurements taken at home can provide your doctor with valuable information on your normal blood pressure values under actual everyday conditions.

The European Society of Hypertension (ESH) has set up the following standard blood pressure values when measured at resting pulse:

Blood pressure (mmHg)	Too low	Normal values	Grade 1 & 2 hypertension	Grade 3 hypertension
SYSTOLIC = systole = upper value	below 100	up to 139	140-179	≥ 180
DIASTOLIC = diastole = lower value	below 60	up to 89	90-109	≥ 110

[†] European Society of Hypertension

Safety information and important precautions


- To ensure accurate measuring results, carefully read the complete instructions for use.
- This product is intended for household use only. Keep product and batteries away from children.
- People suffering from cardiac arrhythmia, vascular constriction, arteriosclerosis in extremities, diabetes or users of cardiac pacemakers should consult their doctor before measuring their blood pressure themselves, since deviations in blood pressure values may occur in such cases.
- If you are under medical treatment or taking any medication, please consult your doctor first.
- The use of this blood pressure monitor is not intended as a substitute for consultation with your doctor.
- Protect from dust and from direct sunlight. Do not open or modify this equipment without authorization of the manufacturer.

Product description

- 1 LCD display
- 2 Power (start/stop) button
- 3 Hose port
- 4 Battery compartment
- 5 Connector
- 6 Arm cuff
- 7 Air hose

Inserting batteries

- Use alkaline batteries, type AA 1.5V as supplied with this product.
- Remove the battery compartment cover and insert four batteries with correct polarity (see symbol in the battery compartment).
- New batteries will provide approximately 400 measurements.

Key rules for accurate blood pressure measurement

- Always take readings at the same time of day, ideally in the morning, under the same conditions.
- Do not measure within 30 minutes after smoking, coffee or tea consumption or any form of exertion. These factors will influence the measurement results.
- Always measure on the same upper arm.
- Take off wrist watch and jewelry before fitting the cuff on the measuring arm.
- Wait approximately 3 minutes before repeating a measurement.
- Do not wrap the cuff over a jacket or sweater. This will affect the measurement.
- Remove tight-fitting clothing from your arm.

Applying the arm cuff


Fig. 1


Fig. 2


Fig. 3

1. Plug the air hose into the connector (Fig. 1).
2. The cuff is located correctly when the velcro is on the outside of the cuff and the metal ring does not touch the skin (Fig. 2).
3. Put your arm through the cuff loop. The bottom of the cuff should be approximately 2-3 cm above elbow. The tube should lie over the brachial artery on the inside of the arm (Fig. 3).


Fig. 4


Fig. 5


Fig. 6

4. Pull the cuff so that the top and bottom edges are tightened around your arm (Fig. 4).

5. When the cuff is positioned properly, press the Velcro firmly against the pile side of the cuff.
6. This cuff is suitable for use if the arrow heads fall within the stripe pattern on the other side of the cuff when the cuff is tightened around your arm (Fig. 5).
7. Sit on a chair and place your arm on the table so that the cuff is at the same level as your heart (Fig. 6).

Taking a measurement

8. Sit on a chair with your feet flat on the floor.
9. Wrap the cuff snugly around your arm (see "Applying the arm cuff" section above, Fig. 3-6).
10. Do not move the unit during measurement, or the proper measurement will not be achieved.
11. Press power button  to begin. The measurement starts automatically. While measurement is in progress, do not move or speak.

Note: If necessary to interrupt a measurement, press the power button  at any time. The device immediately cancels the measurement, lowers cuff pressure and automatically turns off.

12. The cuff will first inflate and then deflate. At the end of the measurement there are 2 final beep sounds and the reading is displayed:

Upper (systolic) value
Lower (diastolic) value
Pulse

13. Press the power button  to turn the product off. Otherwise the product will turn off automatically after 1 minute. If at any time you want to stop the measurement, press the power button .
14. Test results will be saved into the memory automatically (see memory function below).

Memory function


Your ExactFit™ 1 blood pressure monitor will automatically store the latest measurement. Storing is done automatically after the completed measurement.

To recall the reading, press and hold the power button  for at least 2 seconds while the unit is switched off. When all the display elements are shown, release the power button and the screen will show the last reading with an **M** symbol.

Battery charge indicator


Batteries almost discharged

When the batteries are approximately 75% used, a battery warning appears. This icon blinks 3 times first when the device is turned on, then goes to measuring mode. The monitor will continue to measure reliably, but new batteries should be obtained.


Batteries discharged—replacements required

 When the batteries are fully discharged, the battery symbol will flash empty along with 3 short beep sounds as soon as the device is switched on. You cannot take any further measurement and must replace the batteries. Refer to the section called "Inserting batteries".

Storage and cleaning

- Do not expose device to either extreme temperatures, humidity, dust or direct sunlight.
- The cuff contains a sensitive air-tight inflatable bladder. Handle this carefully and avoid all types of straining through twisting or buckling.
- Clean the product with a soft, slightly damp cloth. Stains on the cuff can be carefully removed with a damp cloth. Do not use gasoline, thinners or similar solvents. Do not wash or dry clean the cuff. Do not store the product in a place exposed to direct sunlight or high humidity.
- Do not drop the device or treat it roughly in any way. Avoid strong vibrations.

- NEVER open the device! This will make the manufacturer's warranty invalid and could adversely affect the device's accuracy!

Calibration

This product has been calibrated at the time of manufacture. If used according to the instructions for use, periodic recalibration is not required. If at any time you question the measuring accuracy, please contact our service representative (see last page for contact information).

What to do if ...

Error	Description	Potential cause and remedy
	Batteries are low	Get ready to replace batteries soon with new ones.
	Batteries are fully discharged	Insert new batteries or check that batteries are correctly inserted.
«Err 1»	Signal too weak	The pulse signals on the cuff are too weak. Re-position the cuff and repeat the measurement.*
«Err 2»	Error signal	During the measurement, error signals were detected by the cuff, caused for instance by movement or muscle tension. Repeat the measurement, keeping the arm still.
«Err 3»	No pressure in the cuff	An adequate pressure cannot be generated in the cuff. A leak may have occurred. Check that the cuff is correctly connected and is not too loose. Replace the batteries if necessary. Repeat the measurement.
«Err 5»	Abnormal result	The measuring signals are inaccurate and no result can therefore be displayed. Read through the checklist for performing reliable measurements and then repeat the measurement.*
«H!»	Pulse or cuff pressure too high	The pressure in the cuff is too high (over 300 mmHg) OR the pulse is too high (over 200 beats per minute). Relax for 5 minutes and repeat the measurement.*
«Lo»	Pulse too low	The pulse is too low (less than 40 beats per minute). Repeat the measurement.*

* If this or any other problem occurs repeatedly, please consult your doctor.

Error	Description	Potential cause and remedy
Readings from repeated measurements differ considerably.	Blood pressure is a fluctuating value. For healthy adults, deviations of 10 to 20 mmHg are possible.	Reminder: always use the same arm for measurements!
	Cuff is not properly fitted.	Ensure that the cuff is fitted according to the "Applying the arm cuff" section.
	Readings were not consistently taken at heart level.	For each measurement, make sure that the cuff is at heart level.
	Talking, coughing, laughing, moving etc. when measuring will influence the reading.	While measuring, relax, keep still, do not move or speak.
Blood pressure values measured at your doctor's office differ from your measurements.	Doctor visits can frequently cause anxiety.	Record the daily development of the measured values and consult your doctor.
Display remains blank, or is unusual, when instrument is switched on.	Batteries are not installed properly.	Check batteries for correct polarity.
Device frequently fails to measure blood pressure values, or values measured are too low (or too high).	Cuff may not be properly positioned.	Ensure that the cuff is fitted according to the "Applying the arm cuff" section.

Further information

Blood pressure is subject to fluctuations even in healthy people.

Comparable measurements always require the same conditions (quiet conditions)!

In order to receive market clearance from governmental bodies, this device was subjected to strict clinical tests.

Specifications

Model:	BUA5000
Operation principle:	Oscillometric method
Display:	Liquid crystal display
Cuff pressure display range:	0-299 mmHg
Pulse measurement:	40-199 beats per minute
Blood pressure measurement range:	30-280 mmHg
Laboratory accuracy:	±3 mmHg (cuff pressure) ±5% of reading (pulse rate)
Clinical accuracy:	According to EN 1060-1/-3 and ISO 81060-2 with auscultatory reference: <5 mmHg systematic offset <8 mmHg standard deviation
Inflation:	Automatic
Power supply:	4 batteries, type AA 1.5V
Operating temperature/ atmospheric pressure/ humidity:	10 °C to 40 °C (50 °F to 104 °F) 860 to 1060 hPa (0.849 to 1.046 atm) 15% to 90% relative humidity, non-condensing
Storage temperature/ atmospheric pressure/ humidity:	-20 °C to 55 °C (-1 °F to 131 °F) 860 to 1060 hPa (0.849 to 1.046 atm) 15% to 90% relative humidity, non-condensing
Cuff:	Universal cuff (22-42 cm)
Clinical validation:	According to international protocol 2002 (ESH)
Service life:	5 years
Warranty:	2 years

If device is not used within specified temperature and humidity ranges the technical accuracy cannot be ensured.

Type BF applied parts

Operating temperature

See instructions for use

Storage temperature

Keep dry

Subject to change without notice.

Internally powered equipment.

Continuous operation.

This device conforms to the following standards:
 EN 60601-1: Medical electrical equipment –
 Part 1: General requirements for safety
 EN 1060-1: AMD 1 Non-invasive sphygmomanometers –
 Part 1: General requirements
 EN 1060-3: Non-invasive sphygmomanometers –
 Part 3: Supplementary requirements for electro-mechanical blood pressure measuring systems
 EN 60601-1-2: Medical electrical equipment –
 Part 1-2: General requirements for basic safety and essential performance –
 Collateral standard: Electromagnetic compatibility – Requirements and tests
 EN 60601-1-11: Medical electrical equipment –
 Part 1-11: General requirements for basic safety and essential performance –
 Collateral standard: Requirements for medical electrical equipment and medical electrical systems used in the home healthcare environment

This product conforms to the provisions of the
 EC directive 93/42/EEC (Medical Device Directive).

This product carries the CE mark and is manufactured in conformity with the
 RoHS Directive 2011/65/EU.

MEDICAL ELECTRICAL EQUIPMENT needs special precautions regarding EMC.
 For detailed description of EMC requirements please contact a service
 representative (see Package Insert or consult www.hot-europe.com). Portable
 and mobile RF communications equipment can affect MEDICAL ELECTRICAL
 EQUIPMENT.

This product contains batteries and recyclable electronic waste.
 To protect the environment, do not dispose of it in the household
 waste, but take it to appropriate local collection points.

Guarantee

Consumer Card available on our website at
www.hot-europe.com/after-sales.

For support requests, please visit www.hot-europe.com/after-sales or
 find service contact information on the last page of this manual.

Production date is given by the LOT located in the back of the device.
 The first 3 numeric digits after the LOT No represent the day of the year of
 manufacture. The next 2 numeric digits represent the last two numbers
 of the calendar year of manufacture and the letters at the end designate
 the manufacturer of the product. E.g. LOT 15614ABC: this product was
 made on the day 156, year 2014 at manufacturer with identifier ABC.

Guidance and manufacturer's declaration – electromagnetic emissions

The BUA5000 is intended for use in the electromagnetic environment specified below.
 The customer or the user of the BUA5000 should assure that it is used in such an
 environment.

Emissions test	Compliance	Electromagnetic environment – guidance
RF emissions CISPR 11	Group 1	The BUA5000 uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	Complies
Harmonic emissions IEC 61000-3-2	Not applicable	The BUA5000 is solely battery powered.
Voltage fluctuations/ flicker emissions	Not applicable	

Non-life support equipment separation distance calculation (3Vrms / 3V/m compliance)

Rated maximum output power of transmitter (W)	Separation distance according to frequency of transmitter (m)		
	150kHz to 80MHz in ISM bands $d = \left[\frac{3,5}{V_1} \right] \sqrt{P}$	80MHz to 800MHz $d = \left[\frac{3,5}{E_1} \right] \sqrt{P}$	800MHz to 2.5GHz $d = \left[\frac{7}{E_1} \right] \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.37	0.37	0.74
1	1.17	1.17	2.33
10	3.69	3.69	7.38
100	11.67	11.67	23.33

Guidance and manufacturer's declaration – electromagnetic immunity			
The BUA5000 is intended for use in the electromagnetic environment specified below. The customer or the user of the BUA5000 should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±6kV Contact ±8kV Air	Complies	Floors should be wood, concrete, or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Radiated RF IEC 61000-4-3	3V/m 80MHz to 2.5GHz	Complies	Field strengths outside the shielded location from fixed RF transmitters, as determined by an electromagnetic site survey, should be less than 3V/m.
Conducted RF IEC 61000-4-6	3Vrms 150kHz to 80MHz	Not applicable (no electrical cabling)	Interference may occur in the vicinity of equipment marked with the following symbol:  Separation distance calculation provided above. If a known transmitter is present the specific distance can be calculated using the equations.
Electrical fast transient IEC 61000-4-4	±2kV power line ±1kV I/O lines	Not applicable	The BUA5000 is solely battery powered.
Surge IEC 61000-4-5	±1kV differential ±2kV common	Not applicable	
Power frequency magnetic field IEC 61000-4-8	3A/m	Complies	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	>95% dip 0.5 cycle 60% dip 5 cycles 70% dip 25 cycles 95% dip 5 sec.	Not applicable	The BUA5000 is solely battery powered.

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