- de Testanleitung für die Eigenanwendung des Schnelltests **Preventis SmarTest® Vitamin D** zur quantitativen Bestimmung von 25(OH)-Vitamin D aus Kapillarblut
- en Instructions for carrying out self-testing with the **Preventis SmarTest® Vitamin D** rapid test for the quantitative determination of 25(OH)-vitamin D in capillary blood
- fr Mode d'emploi de l'outil d'autodiagnostic rapide **Preventis SmarTest® Vitamin D** pour l'évaluation du taux de 25(OH)-vitamine D à partir de sang capillaire
- it Istruzioni per l'esecuzione del test rapido **Preventis SmarTest® Vitamin D** volto all'autodeterminazione quantitativa del livello di 25-OH vitamina D da sangue capillare
- tr Kılcal damardaki kandan 25(OH)-Vitamin D'nin niceliksel olarak belirlemesine yönelik, kendi kendinize kullanabileceğiniz **Preventis SmarTest® Vitamin D** hızlı testi için kullanım kılavuzu

Immundiagnostik AG
Stubenwald-Allee 8a
64625 Bensheim, Germany

Vertreib-/Distribution/Distribution/Distribuzione/Dağıtım:
Preventis GmbH
Stubenwald-Allee 8a
64625 Bensheim, Germany

REF PST80100EP04

Temperaturbegrenzung Temperature limitation Limite de température Limite di temperatura Sıcaklık sınırı	Hersteller Manufacturer Fabricator Produttore Üretici	In-vitro-Diagnostikum: nur zur Anwendung außerhalb des Körpers In vitro diagnostic device: only to be used outside of the body Diagnostic in-vitro: utilisation à l'extérieur du corps uniquement Diagnos in-vitro: solo per uso esterno In vitro Tanı: Sadece vücut dışı kullanımı içindir
Artikelnummer Catalogue number Réf. Référence Codice articolo Ürün numarası	Chargennummer Lot number Numero di lotto Número lote Parti numarası	Gebrauchsanweisung beachten Read user instructions Lire la notice Rispettare le istruzioni per l'uso Kullanma talimatına uyun
Nicht zur Wiederverwendung Do not reuse À usage unique Non riutilizzare Tekrar kullanılmı için uygun değildir	Verwendbar bis Expiry date Utiliser avant le Scadenza Son kullanma tarihi	Inhalt ausreichend für «-» Prüfungen Contains sufficient for «-» tests Contenu suffisant pour «-» tests Contenuto sufficiente per «-» test İçerik «-» adet testi içindir
Bitte wenden Please see reverse Voir au verso Si prega di vedere il retro Lütfen arka yüzü bakın	Literatur/References/Publications/Riferimenti bibliografici/Literatür: 1. Bischoff-Ferrari HA (2014) Optimal serum 25-hydroxyvitamin D levels for multiple health outcomes. In Sunlight, Vitamin D and Skin Cancer (Advances in experimental medicine and biology). Reichrath J (ed) pp 500–25. New York: Springer-Verlag. 2. Gupta AK, Jamal V, Sakul & Malhotra P (2014) Hypervitaminosis D and systemic manifestations: A comprehensive review. J. Int. Med. Sci. Acad. 27: 236–237. 3. Kimball SM, Mirhosseini N & Holick MF (2017) Evaluation of vitamin D3 intakes up to 15,000 international units/day and serum 25-hydroxyvitamin D concentrations up to 300 nmol/L on calcium metabolism in a community setting. Dermatolendocrinol. 9: e1300213.	

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Usage
Preventis SmarTest® Vitamin D is a test system for the self-test for the quantitative determination of the quantitative determination of 25(OH)-vitamin D in capillary blood.

Introduction
“Vitamin D” is the term for a group of molecules that perform many different important functions in the body. Vitamin D is either absorbed through food, or is produced in the skin when it is exposed to the sun. In the human body, vitamin D is stored as 25(OH)-vitamin D. 25(OH)-vitamin D accounts for about 95% of the total amount of vitamin D, and is therefore determined as a marker for the supply of vitamin D. A severe vitamin D deficiency is associated with lower bone density, muscle disorders, cardiac and circulatory diseases, disorders of the nervous system and diabetes. A vitamin D deficiency can facilitate the development of cancer, diabetes, hypertension, myocardial infarction and autoimmune diseases. A regular check of vitamin D levels is especially recommended in the case of risk groups such as infants, elderly people, people with darker skin color, and people who spend very little or no time outside.

Test principle
The **Preventis SmarTest® Vitamin D** test system consists of an immunological rapid test that detects 25(OH)-vitamin D via specific antibodies, and a smartphone app that analyzes the result of the rapid test quantitatively based on the color intensity of the control and test lines. The result is displayed on the smartphone used and saved.

Test characteristics
The **Preventis SmarTest® Vitamin D** has a measuring range of 5–100 ng/ml 25(OH)-vitamin D. The test detects 25(OH)-vitamin D₁, 25(OH)-vitamin D₂, 1,25(OH)₂-vitamin D₃ and 1,25(OH)₂-vitamin D₃. There is no cross-reactivity to 1,24(OH)₂-vitamin D₃, 1,24(OH)₂-vitamin D₃ or 3-epi-(OH)-vitamin D₃. The test result of **Preventis SmarTest® Vitamin D** is not affected by bilirubin, triglycerides, cholesterol, vitamin B₁₂ or vitamin C. For repeated measurements with a rapid test rate, the coefficient of variation is 25,8 %, for repeated measurements with different rapid test batches, it is 3–37%. A comparative measurement of 27 patient samples with **Preventis SmarTest® Vitamin D** and the reference method LC-MS/MS (liquid chromatography mass spectrometry/mass spectrometry) showed a very good correlation between the two methods (87,5 %).

Materials
Materials provided

	Test device, individually packed		Camera test card	TEST	TEST	Test instructions
Sampling set:						
	Lancet LANCET		Capillary tube PIP		Buffer tube BUF	
	Cleaning cloth WIPE		Plaster PLASTER			

Additional materials required
Please note: In order to analyze the rapid test, you need the **Preventis SmarTest® Vitamin D** app. This can be downloaded from the Google Play Store or the Apple App Store. The supported smartphones (system requirements) can be found here: <https://vityd.preventis.com>.

Storage and stability
The **Preventis SmarTest® Vitamin D** test system should be stored at 4–30 °C. Do not use the components of the test system after the expiration date. Protect the test from heat and direct sunlight. The rapid test should be used immediately after opening the packaging. Only open the packaging when prompted to do so by the **Preventis SmarTest® Vitamin D** app.

General notes
Please allow around 20 minutes to run the test, without interruption. Ensure that your smartphone has sufficient power (at least 20% battery level). While running the test, do not switch off your smartphone or switch it to silent, as the timer in the app will emit a beep.

- Precautions**
- For *in-vitro* diagnostic use only, i.e. only to be used outside of the body.
 - Do not use the test after the expiry date.
 - Do not use the test device if the packaging is damaged.
 - Do not drink the content of the sample tube, and avoid direct contact with the liquid. In the event of skin contact, immediately wash the affected area thoroughly.
 - The test is designed for use at room temperature. Do not use the test at higher temperatures.
 - The test device must be used immediately after opening the packaging. Only open the packaging when prompted to do so by the **Preventis SmarTest® Vitamin D** app.

Installation of the Preventis SmarTest® Vitamin D app
Install the **Preventis SmarTest® Vitamin D** app from the Google Play Store or from the Apple App Store. The supported smartphones (system requirements) can be found here: <https://vityd.preventis.com>.


- Setting up the Preventis SmarTest® Vitamin D app**
Before initial use, the **Preventis SmarTest® Vitamin D** app has to be set up. To do this, carry out the following:
- Start the **Preventis SmarTest® Vitamin D** app.
 - Read the General Terms and Conditions carefully and confirm that you agree to them by clicking “Continue”.
 - Carry out the camera test. In order to do this, follow the instructions in the app. Carry out the camera test in the same location that you will later use to analyze the rapid test, in order to include the light conditions in the test. **Note:** If the camera test repeatedly fails, check the camera lens for scratches or dirt, and remove them if necessary. Also check the camera card test for possible damage, and ensure that there is even lighting. If the test continues to fail, please contact your **Preventis SmarTest® Vitamin D** supplier.
 - Assign a four-digit personal identification number (PIN) and confirm this by entering it again. You will need your PIN every time you start the app.
 - The **Preventis SmarTest® Vitamin D** app is now set up and ready to use right away.

Functionality of the Preventis SmarTest® Vitamin D app
In the main menu you have the following options:

- “Prepare test” provides information on the preparation of the test, and guides you through the sampling process, running the test, and test analysis, with short illustrated sequences.
- “Results” displays the stored results.
- The “gear symbol” takes you to the settings. Here, amongst other things, you can change the PIN, redo the camera check and reset the app. You can also see detailed instructions here regarding the **Preventis SmarTest® Vitamin D**.


Preparation
Please allow approximately 20 minutes to complete the test. Ensure that the materials provided are at room temperature. Prepare the packaged test device and all other test components. Make sure the tip of the capillary does not come into contact with anything. **Important:** Do not remove the test device from the packaging. Check the internet connection to your smartphone and turn the sound on. Open the **Preventis SmarTest® Vitamin D** app and select “Prepare test”. Then select “To test” to begin the sampling process and completing and analyzing the test.

Sampling
The side regions of the finger tips of your middle or ring finger are particularly suitable for taking a capillary blood sample. Rub the selected finger until it is warm, with good blood circulation. Clean the area with the cleaning cloth, then allow the cleaning liquid to dry completely. Rotate the lancet securing device and remove it. Hold the lancet against the clean finger, and press the yellow button. The lancet is activated, and punctures the finger. Rub the finger and then press it slightly in the direction of the puncture point, until a sufficiently large drop of blood has formed. Hold the capillary horizontal or below the drop of blood and wait until the drop of blood is drawn into the capillary by capillary force alone. Fill the capillary tube completely in this manner with blood without interruption. If air bubbles are visible in the capillary tube, the procedure must be repeated with a new capillary tube. Stick the filled capillary into the buffer tube, piercing the film. Press the capillary tube and buffer tube tightly together until the bridge of the capillary tube completely rests on the edge of the buffer tube (Fig. 1). If the cap comes loose from the capillary when pressing together, please put it back on. Shake the buffer tube and capillary tube vigorously until the blood and buffer have mixed well. You can now cover the puncture point on the finger with a plaster.


Test procedure
Remove the test cartridge from its packaging. Unscrew the cap from the buffer tube/capillary tube. Applying even pressure to the tube, place 3–4 **hanging drops** of the liquid onto the sample application window of the test device (Fig. 2). Start the timer of the **Preventis SmarTest® Vitamin D** app immediately, by pressing “Start timer”. The incubation time of the sample is 15 minutes.

Test analysis
After the incubation period, a beep will be heard. You now have 2 minutes to analyze the rapid test using the app. First check to see whether there is a control line in the results window of the test device. The position of the control line in the results window is indicated by a “C” on the test device. If no control line is visible, the test is not valid and cannot be analyzed! In such an event, the measurement must be repeated, with a new test. The test line is indicated by the letter “T”. The test line may or may not be visible depending on the 25(OH)-vitamin D concentration in the sample (Fig. 3).

When analyzing the rapid test, place the test device on a bright, flat surface. Look for similar light conditions as those present in the camera test. Avoid shadows, strong light from the side and direct sunlight. In the app, press “Scan test”. You will see a template of the test device on your screen. Align the template on your screen with the test cartridge in front of you. While doing this, hold your smartphone horizontally, parallel to the test cartridge and avoid tilting it. Hold it in this position briefly until the **Preventis SmarTest® Vitamin D** app automatically triggers the camera several times, and goes to the analysis screen. After a successful calculation, the **Preventis SmarTest® Vitamin D** app shows the reading on the screen of the smartphone. You can return to the main menu of the app by pressing “Continue”.

Interpretation of the reading^{1,2,3}

0-10 ng/ml	The result indicates a severe deficiency.
11-30 ng/ml	The result indicates a deficiency.
31-100 ng/ml	The result indicates an optimal level.
>100 ng/ml	The result indicates an excessive level.

Please note: Interpretation of the reading may vary from this table, depending on individual health and physical conditions. If you have questions regarding the interpretation of your reading or treatment, please contact your doctor or pharmacist.

Disposal
After use, please dispose of all components of the **Preventis SmarTest® Vitamin D** test kit in the residual waste, including camera test card and test instructions.

Test limitations
The **Preventis SmarTest® Vitamin D** test system detects 25(OH)-vitamin D with a high degree of reliability; however, in certain cases it may produce incorrect results. In the case of questionable results, additional clinical tests should be carried out. As with all diagnostic tests, a definitive diagnosis should not be made on the basis of a single result but only after a doctor has performed a full evaluation of your condition.

Short instructions Preventis SmarTest® Vitamin D

- Clean the fingertip with the cleaning cloth and leave to dry.
- Puncture the skin of the fingertip, using the lancet.
- Fill the capillary tube completely with blood.
- Stick the capillary tube into the buffer tube, push them together completely, and shake vigorously.
- Remove the test device from its packaging.
- Unscrew the cap from the buffer tube. Place 3–4 **hanging drops** of the sample solution into the oval sample application window.
- Observe the 15-minute incubation period.
- After the beep, analyze the result within 2 minutes, using the app.
- The result is displayed on your smartphone.

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